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# TABLE OF CONTENTS

## **Part 1. Technical Passport**

1. Safety rules.....	70
2. Intetnded purpose.....	75
3. Unit description.....	76
4. Application instruction.....	80
5. Package contents.....	89
6. Maintenance and battery replacement procedure.....	90
7. Troubleshooting.....	93
8. Technical information.....	97
9. Manufacturer`s guarantees.....	104
10. Manufacturer`s address.....	106

## **Part 2. Application Instruction**

1. General information.....	108
-----------------------------	-----

2. Indications and contraindications.....	109
3. Treatment conditions.....	113
4. Electrostimulation intensity.....	115
5. Stimulation methods.....	117
6. Indication of modes.....	118
7. Recomendation on selection of treatment areas.....	125

<b>Annex 1. Brief reference atlas on the application of DENS for treatment of common diseases.....</b>	191
--	-----

<b>Annex 2. Accessories used.....</b>	211
---------------------------------------	-----

<i>Warranty Maintenance Form.....</i>	217
---------------------------------------	-----

<i>Acceptance certificate.....</i>	220
------------------------------------	-----

# Thank you for using the DiaDENS-T unit.

In order to make the usage of the unit effective and safe, please read carefully all section of this manual.

# PART 1

EN

## Technical Passport

# 1. SAFETY RULES

**Please read carefully all information marked with this sign.** This information is important for ensuring safety and effective use of the unit.

The information contained in this Operating Manual is important for your safety and for correct application and maintenance of the unit.

The unit is safe in operation as it uses the low voltage power supply isolated from the working part of the unit (type B device with a type F working part).

The unit may not be used to treat patients with implanted electronic devices (such as the heart pacemaker) and to treat patients with idiosyncrasy to electric current.

The unit may not be used in the area of direct anterior projection of the heart.

During stimulation the patient should not be connected to any high frequency electric device. Simultaneous use of the unit and other electric devices may cause burns and possible damage of the unit.

Operation in the vicinity of short-wave or microwave equipment may cause the output parameters of the unit to become unstable.

Only the detachable electrodes manufactured by the supplier may be connected to the unit.

The product contains fragile components. Please avoid shocks.

Avoid long-term exposure of the unit to direct sunrays at a high ambient temperature (over 25 °C). This may cause failure of the indicator on the unit or destroy its batteries.


The unit is not waterproof. Protect it from moisture.


All maintenance works on the unit must be performed by skilled technicians at the manufacturing factory.


**Transportation conditions:** temperature -50 to +50 °C, relative air humidity 30 to 93 %, atmospheric pressure 70 to 106 kPa (525 to 795 mm Hg).

**Storage conditions:** temperature -50 to +40 °C, relative air humidity 30 to 93 %, atmospheric pressure 70 to 106 kPa (525 to 795 mm Hg).

**Operating conditions:** temperature +10 to +35 °C, relative air humidity 30 to 93 %, atmospheric pressure 70 to 106 kPa (525 to 795 mm Hg).

**Warning!** *If the unit was stored at the ambient temperature lower than 10 °C, before operating keep it under normal climatic conditions for at least two hours.*

**Recycling:** the unit packaging materials are not hazardous for the environment; they can be recycled.

Separate collection of electric and electronic equipment.

**Warning!** *The unit contains precious materials that may be used repeatedly after disposal with consideration of the environmental protection regulations. They must be submitted to specially designated collection and recycling facilities (please consult in the relevant services in your locality).*

## 2. INTENDED PURPOSE

The DiaDENS-T unit is used for electrostimulation in a broad range of frequencies with selectable individual treatment programs.

The DiaDENS-T unit is designed for medical and preventive treatment facilities and individual application in accordance with the instructions of the attending physician.


**Compliance:**


*This medical device bears the CE marking in compliance with Directive 93/42/EEC concerning Medical Devices.*

### 3. UNIT DESCRIPTION

EN

Picture 1. Arrangement of the unit. Front



Picture 2. Arrangement of the unit. Back

\* The unit contains a duplicate label with the serial number under the batteries (inside the battery compartment) of the unit. Make sure not to remove it and not to delete the serial number as it results in loss of warranty for the unit.

**\*\* Warning!** The connector is only designed for connection of electrodes manufactured by the supplier. Never connect the unit to the computer as it may cause failure of the computer and expose the user to the risk of electrocution.

Picture 4. Designation of symbols on the unit display

## 4. APPLICATION INSTRUCTIONS

### 4.1. Procedure conditions

No special conditions are required for the procedure. The room should be dry and warm. During the procedure the patient should be positioned comfortably, seated or lying. After the procedure it is recommended to remain relaxed for 10–15 minutes.

### 4.2. Switching on the unit

Press the  button; the unit should turn on and switch to the standby mode.

### 4.3. Mode selection

To select the desired mode, use the **Mode**  and **Mode**  buttons.

The modes in the menu of the unit are located in the following sequence: Screening →

MED → 20 Hz → 60 Hz → 77 Hz → 140 Hz → 200 Hz. When the unit is switched on, the default mode setting is 77 Hz.

#### 4.4. Stimulation power setting

In order to set the desired power, use the **Power**  and **Power**  buttons:

- each time the button is pressed, the power is decreased or increased by one unit;
- for continuous increase or decrease of the stimulation power press and hold down the button. When the desired power level is reached, release the button.

The stimulation power level setting is shown in the left top corner of the display as a number symbol and on the linear horizontal scale. If the detachable electrode is connected, the power is reset to 0.

#### 4.5. Timer

If you need to set a fixed stimulation time, use the Timer function intended for operation in the 200 Hz, 140 Hz, 77 Hz, 60 Hz and 20 Hz modes and is activated by pressing the

button.

By default the Timer function is switched off; in this case the maximum procedure duration is 30 minutes.


Each time you press the  button, the stimulation time is set at 5, 10, 15 and 20 minutes and the Timer Activated indication () is shown on the display.

Next pressure of the button (after a 20-minute interval) switches off the Timer function; the Timer Activated indication on the display disappears.

#### 4.6. Volume setting

The unit can operate in a mute mode or you can select one of the two sound alarm volume settings. The volume level is set using a combination of the , **Power**  and **Power**  buttons.

To select the sound alarm volume, press a combination of  and **Power**  buttons simultaneously to increase the volume or press the combination of  and **Power**  buttons simultaneously if you want to decrease the volume or switch off the sound.

The volume level is indicated on the display of the unit: no picture; ; ; .

#### 4.7. Forced Therapy Function

In the Forced Therapy mode, the unit can operate in the therapy modes irrespective of the presence and the quality of the contact of the electrodes with the skin surface. The mode is primarily used when detachable electrodes are connected, when it is not always possible to control the quality of the contact.

To activate the Forced Therapy function, press simultaneously the **Power**  and  buttons; the mode is deactivated by simultaneously pressing the **Power**  and  buttons.

buttons. When the Forced Therapy function is activated, the contact symbol  on the display starts flashing.

**Warning!** *With the Forced Therapy Function activated, the unit is not switched off automatically when there is no contact between the electrodes and the skin, even if the keyboard buttons are not used. Do not leave the unit with the Forced Therapy function activated as it will cause discharge of the battery before the indicated time.*

#### **4.8. Switching off the unit**

Press the  button and hold it down until the unit switches off.

#### **4.9. Preparation of the unit for the procedure.**

Before using the unit, make sure that there are no impact marks on its case. The liquid crystal display of the unit is not damaged. Make sure that there are no traces of foreign

materials or marks of damage on the electrode.

**Warning!** Do not use the unit if it is damaged.

EN

**Warning!** If the unit was stored in a cold place or a place with high humidity, keep it at room temperature for at least 2 hours.

**4.9.1.** Open the consumer package and take out the unit.

**Warning!** Before the first use, wipe the electrodes with a disinfectant solution (for example, 3% hydrogen peroxide solution) for hygienic purposes.

**4.9.2.** If there are no batteries in the battery compartment, install the batteries.

**4.9.3.** When installing the batteries, make sure that polarity is correct.

## 4.10. Procedure

**4.10.1.** Prepare the unit for the procedure (see paragraph 4.9).

**4.10.2.** Switch on the unit, using the  button.

**4.10.3.** If necessary, set the sound alarm volume as specified above and fixed duration of the procedure, using the Timer function (  ). If required, connect one of the detachable electrodes to the respective plug (picture 1) specified in Annex 2. In case of a small skin contact area or in case of a dry skin turn on the Forced Therapy function on the unit as per paragraph 4.7. Carry out the procedure as specified below. Upon completion of the procedure turn on the Forced Therapy function.

**Warning!** Do not connect to the unit detachable electrodes not provided by the manufacturer.

- 4.10.4.** Set the required stimulation mode, using the **Mode**  and **Mode**  buttons.
- 4.10.5.** To start the procedure place the electrodes on the skin in the stimulation area. The contact symbol  will be displayed.
- 4.10.6.** Gradually increase the power until it reaches the desired level (minimum, comfortable or maximum). The stimulation time will start.

**Warning!** *It is not recommended to carry out the procedure at an intolerable power level!*

If the Timer function deactivated, the procedure time count in the order of increase starts on the digital indication display.

If the Timer function is activated, time count on the digital indication display starts in the order of decrease from the set timer value.

- 4.10.7.** When the maximum procedure duration time (30 minutes) is reached or earlier (if the timer is activated), the alarm sounds and electrostimulation stops. To continue stimulation, either remove the unit from the skin for a short time or change the electrostimulation power or mode. The alarm also sounds when the electrode loses contact with the skin.
- 4.10.8.** After completion of the procedure switch on the unit and clean the electrodes with soft pile-free tissue moistened with a standard disinfectant (such as 3% hydrogen peroxide).

**Warning!** *Inadequate cleaning of the electrodes may lead to allergic reactions or cross-infection in case of use of the unit by several patients.*

**Warning!** *In case of idiosyncrasy, in case of aggravation of the patient's condition during the procedure or skin irritation after the procedure, stop the stimulation and consult a physician.*

## 5. PACKAGE CONTENTS

The package contains the following items:

Name	Number, pieces
DiaDENS-T unit	1
Operation manual combined with the application instructions and technical passport	1
Type LR6/AA battery	2*
Consumer packaging	1

---

\* *The number of batteries may be increased.*

## 6. MAINTENANCE AND BATTERY REPLACEMENT PROCEDURE

**6.1.** Maintenance should include the following operations:

**6.1.1.** External inspection of the unit. Make sure that there are no impact marks on the unit as this may lead to its incorrect operation.


**6.1.2.** Clean the integrated electrode and the detachable electrodes before and after use. The lack of adequate cleaning may cause allergic reactions or skin infection, especially in case of use of the unit by several patients. Use standard disinfectants (such as 3% solution of hydrogen peroxide) and soft pile-free tissues.

**6.2.** If the unit will not be used for a long period of time (more than 3 days), remove the batteries from the battery compartment. This will prevent their early discharge.


- 6.3.** If the battery indicator on the LCD starts flashing, replace the batteries. To do this:
- open the battery compartment as shown in picture 5;
  - remove the used batteries;
  - install new batteries, making sure that the polarity is correct\*;
  - close the battery compartment cover on the unit.

---


\* *Use only the batteries provided for this unit – LR6/AA type with the rated voltage of 1.5V.*

1. Remove the battery compartment cover

2. Install new batteries, making sure that the polarity is correct

3. Close the battery compartment cover

## 7. TROUBLESHOOTING

The table below describes the situations that may be interpreted as fault and can be remedied without external help. In case of other defects contact the manufacturer's representative (see Section 9 "Manufacturer's Guarantees"), do not try to remedy them.

*Possible faults and remedies*

<b>Problem</b>	<b>Possible cause</b>	<b>Remedy</b>
The unit does not turn on when the  button is pressed.	No batteries	Install new batteries (see Maintenance and Battery Replacement Procedure)
The unit does not turn on when the  button is pressed.	Battery voltage less than 2.1 V	Replace batteries (see Maintenance and Battery Replacement Procedure)
When the unit is turned on, alarm sounds and the unit switches off automatically	Battery voltage less than 2.2 V	Replace batteries (see Maintenance and Battery Replacement Procedure)

Problem	Possible cause	Remedy
When the unit is turned on, alarm sounds and the unit switches off automatically	Time interval between de-energizing and consequent energizing of the unit less than 3 sec.	Keep the unit de-energized for at least 3 sec. and then turn it on again
The unit switches off automatically or the battery symbol flashes and intermittent alarm sounds	Battery voltage less than 2.1 V	Replace batteries (see Battery Replacement Procedure)
The unit does not switch to the stimulation mode when the integrated electrodes are in contact with the skin surface; no skin contact symbol on the display	Zero power setting	Increase power setting
	Contact area between the integrated electrodes and the skin not sufficient	Apply the integrated electrodes of the unit to the skin unit, making sure the contact is sufficient. If necessary, moisten the contact area or turn on the Forced Therapy function, using a combination of the <b>Power</b>  and  buttons.

Problem	Possible cause	Remedy
No sound alarm	The sound is off	Turn on the sound, using the  , and <b>Power</b>  buttons
Quick consumption of batteries	Poor quality batteries	Use good quality batteries (alkaline batteries are recommended) or rechargeable batteries of the corresponding type and size with the voltage not exceeding 1.5 V
	Forced Therapy function activated	Deactivate the Forced Therapy function, using a combination of the <b>Power</b>  and  buttons.

Problem	Possible cause	Remedy
The unit does not switch off automatically when there is no contact between the electrodes and the skin surface and when its controls have not been used for more than 2 minutes	A detachable electrode is connected to the unit.	Turn off the unused detachable electrode
	Forced Therapy function activated	Deactivate the Forced Therapy function, using a combination of the <b>Power</b>  and  buttons.

## 8. TECHNICAL INFORMATION

**8.1.** Maximum consumption current of the unit at the supply voltage of  $3,0 \pm 0,1$  V;  
— in the stimulation mode — max. 300 mA;  
— in a de-energized condition — max. 10  $\mu$ A.

**8.2.** Weight without batteries — max 300 g.

**8.3.** Dimensions — max. 135×51×46 mm.

**8.4.** The unit turns off automatically if not used for more than 3 min. (except in the Forced Therapy mode).

**8.5.** Operating parameters:



**8.5.1.** Stimulation impulse frequency in the Therapy mode:  $20 \pm 2$  Hz;  $60 \pm 2$  Hz;  
 $77 \pm 2$  Hz;  $200 \pm 5$  Hz;





**8.5.2.** Stimulation impulse frequency in the Screening and MED modes —  $10 \pm 2$  Hz.


### 8.5.3. Impulse parameters at minimum power with the load of $(20 \pm 10\%)$ kOhm:


- amplitude in the positive part – max 40 V;
- amplitude in the negative part – max 10 V.

### 8.6. Relation between the impulse waveform and the load resistance at the medium and maximum power level:

Load resistance	Power level: Half		Power level: full	
	Waveform	$V_{p-p}$	Waveform	$V_{p-p}$
Free output		$\approx 150 \text{ V} \pm 20 \%$		$\approx 150 \text{ V} \pm 20 \%$

200 Ohm		$\approx 15 \text{ V} \pm 20 \%$		$\approx 20 \text{ V} \pm 20 \%$
500 Ohm		$\approx 30 \pm 20 \%$ $E_{\text{pulse}} \approx 155 \text{ } \mu\text{J}$		$\approx 40 \text{ V} \pm 20 \%$ $E_{\text{pulse}} \approx 299 \text{ } \mu\text{J}$

### 8.7. Effective current – $I_{\text{eff}}$ for load resistance 500 Ohm:

Mode	Power level: Half	Power level: Full
20 Hz	$I_{\text{eff}}=2,5\text{ mA}$	$I_{\text{eff}}=3,5\text{ mA}$
60 Hz	$I_{\text{eff}}=4,5\text{ mA}$	$I_{\text{eff}}=6,0\text{ mA}$
77 Hz	$I_{\text{eff}}=5,0\text{ mA}$	$I_{\text{eff}}=6,8\text{ mA}$
140 Hz	$I_{\text{eff}}=6.6\text{ mA}$	$I_{\text{eff}}=9.2\text{ mA}$
200 Hz	$I_{\text{eff}}=7,9\text{ mA}$	$I_{\text{eff}}=11,0\text{ mA}$

### 8.8. Electromagnetic radiation

Test	Conformity	Conditions of use
HF radiation CISPR 11	Class B	The DiaDENS-T unit can be used in all kinds of facilities, including application at home

### 8.9. HF radiation resistance

Test	IEC 60601-1-2, testing conditions	Acceptable level
IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz	3 Vrms
IEC 61000-4-3	3 V/m, 80 MHz to 2,5 GHz	3 V/m

### 8.10. Electromagnetic field resistance

Test	Test level	Conformity level	Conditions of use
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 6$ kV contact $\pm 8$ kV air	$\pm 4$ kV contact $\pm 8$ kV air	The floor should be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic materials, relative air humidity should be at least 40 %

Magnetic fields IEC 61000-4-8	3 A/m	3 A/m	The magnetic field parameters should not exceed the limits for commercial buildings and medical facilities
----------------------------------	-------	-------	--

**8.11.** Recommendations concerning determination of the required distance between the DiaDENS-T stimulation unit and radio frequency radiating equipment.

Declared maximum output power of the transmitter P (W)	Radiation frequency and the d distance determination formula (m)		
	150 kHz — 80 MHz $d = 1,2\sqrt{P}$	150 kHz — 800 MHz $d = 1,2\sqrt{P}$	800 MHz — 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

## 9. MANUFACTURER'S GUARANTEES

**9.1.** The manufacturer guarantees failure-free operation of the unit if the operation, transportation and storage conditions are followed.

**9.2.** Lifetime of the unit — 5 years.

If the operating instructions are followed, the lifetime may be much longer than the officially established one.

**9.3.** Warranty operation term — 24 months of the date of sale.

**9.4.** The Seller (Manufacturer) or the organization acting on behalf of the Seller (Manufacturer) under the contract with him shall not be liable for defects if they have emerged after the transfer of the unit to the user due to:

- 1) failure of the user to observe the transportation, storage, maintenance and operation rules provided for by this Manual;

- 2) mechanical damage;
- 3) third parties' actions;
- 4) circumstances of force majeure.

**9.5.** The warranty does not cover the units with destroyed factory seals.

**9.6.** In case of failure of the unit or its defect during the effective term of the warranty and if the owner discovers that certain components are missing, the owner of the unit shall send to the manufacturer or his representative the unit and a repair (replacement) request stating his full name, address and telephone number and containing a short description of the defect and the conditions in which and the date when it was discovered.

## 10. MANUFACTURER'S ADDRESS:

RC ART LLC

Academician Postovsky St. 15, Ekaterinburg 620146, Russia

Tel.: +7 (343) 267-23-30

<http://www.denascorp.ru>

e-mail: [corp@denascorp.ru](mailto:corp@denascorp.ru)

***Official Representative in the European Union:***

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## Application Instructions

# 1. GENERAL INFORMATION

Stimulation of reflexogenous zones and points with impulse currents for preventive treatment, therapy and recovery of bodily functions is an independent method of treatment and can be applied both as part of comprehensive treatment and as individual therapy.

Numerous studies show that electroneurostimulation (DENS) is based on multi-level reflectory and neurochemical reactions that trigger a cascade of regulatory and adaptive mechanisms of the body, which leads to elimination of pain syndromes and blood circulation improvement, has an anti-inflammatory effect, activates production of biologically active substances and metabolism in the tissues and normalizes muscular and vascular tonus. Dynamic electroneurostimulation promotes general health improvement and fitness.

## 2. INDICATIONS AND CONTRAINDICATIONS

The DiaDENS unit may be applied for treatment of patients any age from newly born children to elderly people.

### ***Indications:***

- acute and chronic pain syndromes;
- injuries;
- treatment of functional disorders in case of diseases of the respiratory apparatus, ENT organs, the cardiovascular, musculoskeletal, nervous, endocrine and urinogenital systems, the eyes and the skin;
- rehabilitation after diseases, surgery and injuries;
- enhancement of adaptability of the body in case of exposure to adverse factors and in the conditions of hard physical and mental work, in case of physical and intellectual overstress, chronic fatigue syndrome, difficulties waking up in the morning and sleepiness during the day, in case of dyskoimesis in the evening

and insomnia, increased irritability and depression and for prevention of catarrhal conditions.

### ***Treatment effect***

When using the unit, one should remember that the most apparent therapeutic effect is achieved after a course of treatment. Its duration is 6 to 8 procedures for some diseases, 8 to 12 procedures for other diseases or in some cases — 14 to 20 procedures. The combined therapeutic effect of individual procedures also creates an after-effect of the treatment course that continues after its completion and the long-term results are in some cases even more favourable than the direct ones. But one should not forget that long-term continuous exposure to one physical factor may cause the body to adapt to it, which may considerably decrease the treatment effect.

Moreover, the stochastic nature of physiological processes leads to the inadequate effect of or the lack of favourable response to therapeutic physical factors in some patients (1–5%).

**Warning!** Application of the unit does not rule out the use of other treatment methods indicated by the attending physician.

**Warning!** In case of lack of positive changes or aggravation of the patient's condition consult the attending physician.

### **Contraindications**

#### *Absolute:*

- idiosyncrasy to electric current;
- an implanted heart pacesetter;

#### *Relative\*:*

- status epilepticus;
- neoplasia of any etiology and localization;

---

\* In cases listed above consent of the physician is needed before indicating the treatment.

- acute fever of unknown etiology;
- venous thrombosis;
- acute mental, alcoholic or narcotic excitement

**Warning!** *If you have any of the above listed contraindications, consult your physician before using this device.*

**Warning!** *In case of idiosyncrasy, if you feel worse or the pain increases during the procedure, stop the stimulation and consult a physician immediately.*

### 3. TREATMENT CONDITIONS

**Warning!** Sudden pain of any localization may be the first and frequently the only sign of a serious disease. That is why if pain attacks happen for the first time and repeat and the pain becomes more intense, consult the physician as soon as possible.

DENS procedures do not require any special conditions. The procedures can be performed without attendance of a physician or a nurse or by an operator who will carry out the treatment in the zones and points that the patient cannot reach without help. During the procedure, the patient may be seated in an armchair or lie in a comfortable position. It is recommended to take a 10-15 minute rest after the procedure.

These instructions describes the stimulation schemes used in case of the most common diseases. Additional information on possible applications of DENS can

be found in the Manual on Dynamic Electroneurostimulation Using DiaDens-T and DiaDENS-DT Units\*.

**Warning!** After each procedure clean the electrodes with a standard disinfectant solution (such as 3% hydrogen peroxide solution). The unit should be stored with dry electrodes.

---

\* *Manual on dynamic electrostimulation using DiaDENS and DiaDENS-DT units // Edited by V.V.Chernyshev. – Ekaterinburg, 2005.*

## 4. ELECTROSTIMULATION INTENSITY

Electrostimulation intensity (power) is set individually depending on the patient's individual sensations. For convenience electrostimulation power is divided into three levels (power ranges).

**The first, or minimum, level** — you do not have any subjective sensations or feel a slight vibration in the area under the electrodes. This level is used when intensive stimulation is not recommended, that is for patients with low intensity pain. In case of cardiovascular diseases (arterial hypertension, arterial hypotension and vegeto-vascular dystonia syndrome) the procedure should also start with the stimulation at the minimum power level.

**The second, or comfortable, level** — the patient feels vibration, pleasant tingling or a slight “burning” but without any painful sensations. This level is used for treatment of patients complaining of slight and medium intensity pain. This is the most frequently used power level.

**The third, or maximum level** — the patient has a painful tingling or burning sensation. Stimulation of such intensity may be accompanied by involuntary contraction of the muscles located close to the electrodes.

This level is used for the treatment of a strongly manifested pain syndrome.

**Warning!** *At different stages of treatment the electrostimulation power level can be increased and decreased as the patient's individual tolerance changes and as the pain decreases.*

**Warning!** *The power level should be adjusted depending on the patient's sensations. Do not exceed the pain threshold.*

## 5. STIMULATION METHODS

During the procedure, you can use three stimulation methods — the *stable*, the *labile* and the *combined labile and stable method*.

*The stable method*, (fixed position of the electrodes) is applied for treatment of small areas and areas with skin changes (post-surgery and post-burn scars, oedemas etc.).

When using the *labile treatment method*, move the built-in electrodes smoothly around the treated area so that they remain in constant contact with the body surface, at the speed of 0.5 to 3 cm/sec. You may make linear, spiral, circular or other movements, depending on the area and relief of the treated zone, and only on undamaged skin.

To apply the *combined labile and stable treatment method*, move the electrodes across the skin surface, keeping the built-in electrodes in a fixed position in some points for some time (for example, in maximum pain areas).

## 6. INDICATION OF MODES

### 6.1. Therapy mode **Th**

Stimulation in the Therapy mode is used in case of local pain syndrome, functional disorders and as an emergency care method. Duration of the procedure in the Therapy mode in the area of direct projection of pain and functional disorders is selected based on the following response:

- the complaint is fully eliminated;
- the patient feels considerably better;
- bright redness of the skin or a feeling of warmth in the area under the electrodes.

The following frequencies are available in the Therapy mode:

a) *therapeutic frequency of 140 Hz, 200 Hz* is used in the direct projection area of the complaint. The effect is achieved within the first few minutes and continues from several minutes to 1 hour. To enhance the therapeutic effect, the treatment may be continued at a lower frequency.

Indications: acute pains related to diseases and affection of the musculoskeletal system and the peripheral nervous system damage:

b) *therapeutic frequency 60 and 77 Hz*, is used in the area of direct projection of the complaint and in segmental areas. The effect is achieved within 15–30 minutes and continues up to 1–1.5 hours and more.

Indications: functional disorders and pain syndromes in case of internal diseases and during the post-surgery period;

c) *therapeutic frequency 20 Hz*, is used in the area of direct projection of the complaint, in segmental and universal areas. The maximum effect is achieved within 20 – 60 minutes and continues for several hours.

Indications: functional disorders and pain syndromes in case of internal diseases and diseases of the musculoskeletal system and during the post-surgery period;

## 6.2. MED mode (preventive treatment) **M**


Stimulation in the MED (**M**inimal **E**ffective **D**oze) mode is used in cases of expected hard physical and mental work, in case of physical or mental overstress, chronic fatigue syndrome, difficulties waking up in the morning, sleepiness during the day, inability to concentrate and focus attention, for non-specific preventive treatment of catarrhal diseases during epidemics.

Treatment is carried out in the hand area (he-gu zone) as shown in picture 6. The procedure should be carried out at the comfortable power level (slight tingling). If you do not feel anything, press the **Power** **+** button several times until you feel slight tingling. If you have unpleasant (painful) sensations under the electrode, press the **Power** **-** button several times until you feel slight tingling.

The MED mode can be used one time during each procedure. It is recommended to deliver a course of 10 to 14 procedures.

The unit may be moved to the next area after the mode time is over (sound alarm or indication on the display).

### 6.3. Screening mode **S**

In this mode the user can select the optimum treatment areas for enhancement of the systemic effect of DENS. In this mode the unit detects the location of the latent trigger zones (hidden problem zones) by means of assessment of skin electric resistance growth in the areas adjacent to the selected area.

Procedure: place the electrodes of the unit in the selected area (direct projection of the complaint).

**Warning!** *During operation of the unit in the Screening mode, the electrodes should be placed on the skin in the stable mode, that is you cannot move the electrodes while the unit is operating in the Screening mode.*

When the unit detects contact between the electrodes and the skin surface, it begins the 5 seconds interval, during which is required to detect the skin electric resistance variation in response to the impulse transmitted by the unit. After 5 seconds, the unit

displays the measurement in the form of the index (within the range between 0 and 100 units), for example  $\Delta LT = 8$ . Record the value obtained.

Place the unit in the next skin area within the selected zone. Repeat the manipulation on the whole area.

The skin areas with the maximum  $\Delta LT$  (significantly different from the other numbers) will be the trigger zones.

Stimulate the detected trigger zones additionally in the Therapy mode for 10 – 15 minutes at the frequency of 60 or 77 Hz.

## 7. RECOMMENDATIONS ON SELECTION OF TREATMENT AREAS

**Warning!** *Only the skin can be exposed to stimulation by means of the electrodes.*

Direct stimulation of the lesion focus – in the lumbal area in case of lumbodynia or on the skin in the area around the affected joints in case of pain in the knee joints – is one of the key, simple and rather effective methods of treatment area selection.

To enhance the treatment effect, stimulation of the lesion focus can be combined with stimulation of the respective segmental reflexogenous areas (Annex 1).

It is recommended not to treat more than 2 or 3 zones during one procedure.

Recommended average duration of one procedure:

- for children aged 0 to 12 months – 5-10 minutes;
- for children aged 12 to 36 months – 10-15 minutes;
- for children aged 3 to 5 years – 15-20 minutes;
- for children aged 5 to 12 years – 20-25 minutes;
- for children older than 12 years of age and adults – up to 30 minutes.

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DE

## **ANNEX 1**

Brief reference atlas on the application of DENS  
for treatment of common diseases

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
**ACUTE RESPIRATORY AND ENT  
DISEASES**

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**Nose and the area around it:** 77 Hz — 10 minutes

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
**Area on the posterior surface of the back around the neck**  
20 Hz — 10-15 minutes

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
**TRACHEITIS,  
BRONCHITIS AND  
BRONCHIAL ASTHMA**

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


**Area on the posterior surface of the back around the neck (in case of diseases of the arm joints):** 20 Hz — 7–10 minutes.

**Lumbosacral area (in case of diseases of the leg joints):** 20 Hz — 7–10 minutes

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**Direct projection of the affected joint (marked in red):** 77 Hz (in case of strong pain – 140 or 200 Hz) — until the pain subsides.

**Symmetrical zone (marked in green)** 20 Hz — 10–15 minutes.

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

**Direct projection of pain in the neck or in the back** 77 Hz (in case of strong pain – 140 or 200 Hz) — until the pain subsides.

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**Area on the posterior surface of the back around the neck:** 20 Hz — 7–10 minutes.

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

**Area around the second cervical  
vertebra:** 77 Hz — 5–7 minutes;  
20 Hz — 5–10 minutes

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**Area on the posterior surface of the back around the  
neck:** 20 Hz — 7–10 minutes.

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

**Projection of the affected  
tooth:** 200 and 77 Hz until  
the condition improves

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**Area on the posterior surface of the back  
around the neck:** 20 Hz — 10-15 minutes.

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



**Area around the wound (burn):** 140 and 77 Hz — 15 minutes.

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**Area on the posterior surface of the back around the neck (in case of injury of the upper extremities), lumbosacral area (in case of injury of the lower extremities):** 20 Hz —10 – 15 minutes

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## CONSTIPATION AND ABDOMINAL SWELLING

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

**Projection of the liver and the  
intestines:** 77 Hz — 15 minutes

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

*A set of detachable zonal electrodes **DENS-applicator**. Designed for non-invasive therapy of the pain area, lesion foci and reflexogenous areas (with undamaged skin).*

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**Therapeutic detachable massage electrodes.** Designed for stimulation in the area of the pain, lesion foci, reflexogenous areas and acupuncture points of the body.

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
*Detachable paraorbital electrode **DENAS-goggles**. Designed for preventive therapy and treatment of biologically active points around the eyes in the paraorbital area using the dynamic electroneurostimulation method.*

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**Detachable therapeutic electrode.** Designed for therapy of biologically active points using the dynamic electroneurostimulation method.

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